


Anti-Cdkn2C/P18-Ink4C Polyclonal Antibody For Wb, Elisa - P42773

Item Number: CM0016173



Introduction

Rabbit polyclonal antibody targeting CDKN2C/p18-INK4C, validated for Western blotting and ELISA. This antibody recognizes the cyclin-dependent kinase inhibitor, which regulates cell cycle progression via CDK6/CDK4 interaction.

[Learn More](#)

Parameter	Value
Product Name	CDKN2C/p18-INK4C Rabbit pAb
Remarks/Alias	p18; INK4C; p18-INK4C; CDKN2C/p18-INK4C
Species	Human
GeneID (Human)	1031
GeneID	1031
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-168 of human CDKN2C/p18-INK4C (P42773)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P42773
Protein Weight	18kDa
Shipping	Ice bag