


Anti-Cdkn2C/P18-Ink4C Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P42773

Item Number: CM0020641



Introduction

High-quality rabbit polyclonal antibody against CDKN2C/p18-INK4C, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat samples. Ideal for cell cycle research and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	CDKN2C/p18-INK4C Rabbit pAb
Remarks/Alias	p18; INK4C; p18-INK4C; CDKN2C/p18-INK4C
Species	Human
GeneID (Human)	1031
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-168 of human CDKN2C/p18-INK4C (NP_523240.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42773
Protein Weight	18kDa
Shipping	Ice bag