


Anti-Cdkn2C/P18-Ink4C Rabbit Monoclonal Antibody For Wb - P42773

Item Number: CM0011179



Introduction

Rabbit monoclonal antibody against CDKN2C/p18-INK4C, validated for WB, IF/ICC, ELISA. Recognizes human, mouse, and rat targets. Ideal for cell cycle and cancer research. Recombinant immunogen, 18 kDa.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] CDKN2C/p18-INK4C Rabbit mAb
Remarks/Alias	p18; INK4C; p18-INK4C; 4C
Species	Human
Gene ID (Human)	1031
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human CDKN2C/p18-INK4C (P42773).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42773
Protein Weight	18kDa
Shipping	Ice bag