


Anti-Cdkn2D Polyclonal Antibody For Wb, If/Icc, Elisa - P55273

Item Number: CM0028936



Introduction

High-quality rabbit polyclonal antibody against human CDKN2D (p19-INK4d) for WB, IF/ICC, and ELISA. Reacts with human and mouse CDKN2D. Ideal for cell cycle and tumor suppressor research.

[Learn More](#)

Parameter	Value
Product Name	CDKN2D Rabbit pAb
Remarks/Alias	p19; INK4D; p19-INK4D; CDKN2D
Species	Human
GeneID (human)	1032
GeneID	1032
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-166 of human CDKN2D (NP_001791.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P55273
Protein Weight	18kDa
Shipping	Ice bag