


Anti-Cdkn3 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16667

Item Number: CM0019319



Introduction

High-quality rabbit polyclonal antibody against human CDKN3, validated for Western blot, IHC-P, IF/ICC, and ELISA. Cross-reacts with Human and Mouse. Ideal for cell cycle regulation studies.

[Learn More](#)

Parameter	Value
Product Name	CDKN3 Rabbit pAb
Remarks/Alias	KAP; CDI1; CIP2; KAP1; CDKN3
Species	Human
Gene ID (Human)	1033
Gene ID	1033
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q16667
Protein Weight	24kDa
Shipping	Ice bag