

## Anti-Cdsn Rabbit Polyclonal Antibody For Wb, Elisa - Q15517

Item Number: CM0015097

### Introduction

High-quality anti-CDSN rabbit polyclonal antibody validated for Western blot and ELISA. Shows cross-reactivity with human, mouse, and rat. For epidermal barrier and skin research.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CDSN Rabbit pAb                                                                                                           |
| Remarks/Alias    | PSS; HTSS; PSS1; HTSS1; HYPT2; CDSN                                                                                       |
| Species          | Human                                                                                                                     |
| GeneID (Human)   | 1041                                                                                                                      |
| GeneID           | 1041                                                                                                                      |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human CDSN (NP_001255.3). |
| Source           | Rabbit                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                     |
| Application      | WB, ELISA                                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                         |
| SWISS            | Q15517                                                                                                                    |
| Protein Weight   | 52kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |