

Anti-Cdx2 Rabbit Monoclonal Antibody For Wb And Elisa - Q99626

Item Number: CM0003450



Introduction

High-performance rabbit monoclonal antibody targeting human CDX2 (Q99626). Validated for Western blot and ELISA; cross-reacts with mouse CDX2. Ideal for intestinal epithelium studies, colorectal cancer research, and transcriptional regulation analysis.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] CDX2 Rabbit mAb
Remarks/Alias	CDX3; CDX-3; CDX2/AS; X2
Species	Human
Gene ID (Human)	1045
Gene ID	1045
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human CDX2 (Q99626)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99626
Protein Weight	34kDa
Shipping	Ice bag