

Anti-Cdx2 Rabbit Monoclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - Q99626

Item Number: CM0002453



Introduction

High-quality anti-CDX2 rabbit monoclonal antibody validated for WB, IP, IF-P, IHC-P, and ELISA. Recognizes human, mouse, and rat CDX2, a key homeobox transcription factor in intestinal development and differentiation. Ideal for gut biology and oncology studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] CDX2 Rabbit mAb
Remarks/Alias	CDX3; CDX-3; CDX2/AS; X2
Species	Human
Gene ID (Human)	1045
Gene ID	1045
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CDX2 (Q99626).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99626
Protein Weight	34kDa
Shipping	Ice bag