


Anti-Cdy1B Polyclonal Antibody For Wb, Elisa - Q9Y6F8

Item Number: CM0016053



Introduction

Rabbit polyclonal antibody targeting human CDY1B, suitable for Western blot and ELISA. Detects testis-specific chromatin modifier with histone acetyltransferase activity; cross-reacts with mouse and rat. Ideal for spermatogenesis and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	CDY1B Rabbit pAb
Remarks/Alias	CDY; CDY1B
Species	Human
GeneID	253175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-320 of human CDY1B (NP_001003895.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Y6F8
Protein Weight	60kDa
Shipping	Ice bag