


Anti-Cdyl2 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N8U2

Item Number: CM0031091



Introduction

High-quality rabbit polyclonal antibody specific for Chromodomain Y-like protein 2 (CDYL2/PCCP1). Suitable for Western blot and ELISA, cross-reacts with human, mouse, and rat. Supports nuclear protein and epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	CDYL2 Rabbit pAb
Remarks/Alias	PCCP1; CDYL2
Species	Human
GeneID (Human)	124359
GeneID	124359
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-230 of human CDYL2 (NP_689555.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N8U2
Protein Weight	57kDa
Shipping	Ice bag