

Anti-Ceacam5 (Cea/Cd66E) Rabbit Polyclonal Antibody For Wb, Elisa - P06731

Item Number: CM0018146



Introduction

High-quality CEACAM5 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human carcinoembryonic antigen (CEA/CD66e), ~77 kDa. Ideal for colon cancer research and diagnostic development.

[Learn More](#)

Parameter	Specification
Product Name	CEACAM5 Rabbit pAb
Remarks/Alias	CEA; CD66e; CEACAM5
Species	Human
GeneID (human)	1048
GeneID	1048
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-685 of human CEACAM5 (NP_001278413.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P06731
Protein Weight	77 kDa
Shipping	Ice bag