

Anti-Ceacam5 Monoclonal Antibody For Wb, Ihc-P, Elisa - P06731

Item Number: CM0002657



Introduction

High-quality mouse anti-human CEACAM5 monoclonal antibody validated for WB, IHC-P, ELISA. Recombinant immunogen (aa 35-685). Recognizes CEA/CD66e. Suitable for cancer research, cell adhesion studies, and assay development.

[Learn More](#)

Parameter	Value
Product Name	CEACAM5 Mouse mAb
Remarks/Alias	CEA; CD66e; CEACAM5
Species	Human
Gene ID (Human)	1048
Gene ID	1048
Immunogen	Recombinant fusion protein containing amino acids 35-685 of human CEACAM5
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P06731
Protein Weight	77kDa
Shipping	Ice bag