

Anti-Ceacam5 Polyclonal Antibody For Wb, Elisa - P06731

Item Number: CM0017060



Introduction

Rabbit polyclonal antibody targeting human CEACAM5 (CEA, CD66e). Validated for WB and ELISA; cross-reacts with mouse. Suitable for studying cell adhesion, tumor progression, and colorectal carcinoma research.

[Learn More](#)

Parameter	Value
Product Name	CEACAM5 Rabbit pAb
Remarks/Alias	CEA; CD66e; CEACAM5
Species	Human
Gene ID (Human)	1048
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-702 of human CEACAM5 (P06731)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06731
Protein Weight	77kDa
Shipping	Ice bag