

Anti-Ceacam5 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P06731

Item Number: CM0005188



Introduction

Rabbit monoclonal antibody against human CEACAM5 (CEA/CD66e), validated for Western blot, IHC paraffin, and ELISA. Recombinant immunogen spanning the extracellular domain; ideal for studying cell adhesion, tumor progression, and colorectal carcinoma. Reacts with human.

[Learn More](#)

Parameter	Value
Product Name	CEACAM5 Rabbit mAb
Remarks/Alias	CEA; CD66e; CEACAM5
Species	Human
GeneID (human)	1048
GeneID	1048
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-685 of human CEACAM5 (NP_004354.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P06731
Protein Weight	77kDa
Shipping	Ice bag