

Anti-Ceacam5/Cd66E Rabbit Monoclonal Antibody () For Flow Cytometry - P06731

Item Number: CM0012599



Introduction

-conjugated rabbit anti-human CEACAM5/CD66e (CEA) monoclonal antibody for flow cytometry. Detects a 76 kDa glycoprotein involved in cell adhesion and tumor progression.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CEACAM5/CD66e mAb
Remarks/Alias	CEACAM5; CD66e; CEA; carcinoembryonic antigen-related cell adhesion molecule 5
Species	Human
Gene ID (Human)	1048
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 35-685 of human CEACAM5/CD66e(NP_004354.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06731
Protein Weight	76kDa
Shipping	Ice bag