


Anti-Ceacam6 Monoclonal Antibody For Wb, Elisa - P40199

Item Number: CM0003990



Introduction

Rabbit monoclonal antibody against human CEACAM6 (CD66c) for WB and ELISA. Recombinant immunogen (aa 35-320) targets this 37kDa cell adhesion protein implicated in tumor progression and neutrophil interactions. Swiss-Prot: P40199.

[Learn More](#)

Parameter	Value
Product Name	CEACAM6 Rabbit mAb
Remarks/Alias	NCA; CEAL; CD66c; CEACAM6
Species	Human
Gene ID	4680 (Human)
Immunogen	Recombinant fusion protein containing amino acids 35-320 of human CEACAM6 (NP_002474.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P40199
Protein Weight	37kDa
Shipping	Ice bag