


Anti-Ceacam6 Monoclonal Antibody For Wb, If-P, Elisa - P40199

Item Number: CM0011333



Introduction

Rabbit monoclonal antibody against human CEACAM6 (CD66c), validated for WB, IF-P, ELISA. Cross-reacts with mouse and rat. Suitable for cell adhesion, tumor progression, and neutrophil studies.

[Learn More](#)

Parameter	Value
Product Name	CEACAM6 Rabbit mAb
Remarks/Alias	NCA; CEAL; CD66c; CEACAM6
Species	Human
GeneID (Human)	4680
GeneID	4680
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human CEACAM6 (P40199)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40199
Protein Weight	37kDa
Shipping	Ice bag