


Anti-Ceacam6 Polyclonal Antibody For Wb, Elisa - P40199

Item Number: CM0026644



Introduction

Rabbit polyclonal antibody against human CEACAM6, validated for WB and ELISA. Recognizes human and mouse CEACAM6; suitable for cell adhesion and tumor progression studies.

[Learn More](#)

Parameter	Value
Product Name	CEACAM6 Rabbit pAb
Remarks/Alias	NCA; CEAL; CD66c; CEACAM6
Species	Human
GeneID (Human)	4680
GeneID	4680
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-320 of human CEACAM6 (NP_002474.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P40199
Protein Weight	37kDa
Shipping	Ice bag