

Anti-Cebp Delta/Cebpd Monoclonal Antibody For Wb, Elisa - P49716

Item Number: CM0002205

Introduction

Rabbit monoclonal CEBPD antibody validated for Western blot and ELISA. Cross-reacts with human and mouse. Suitable for immune and inflammatory gene regulation studies. Predicted MW 28 kDa.

[Learn More](#)

Parameter	Value
Product Name	CEBP Delta/CEBPD Rabbit mAb
Remarks/Alias	CELF; CRP3; C/EBP-delta; NF-IL6-beta; CEBP Delta/CEBPD
Species	Human
Gene ID (Human)	1052
Gene ID	1052
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49716
Protein Weight	28kDa
Shipping	Ice bag