

Anti-Cebp Delta/Cebpd Polyclonal Antibody For Wb, Elisa - P49716

Item Number: CM0018269

Introduction

Rabbit polyclonal antibody targeting human CEBP Delta/CEBPD protein, suitable for WB and ELISA applications. Cross-reacts with human, mouse, and rat. Immunogen: synthetic peptide (aa 50-150). Ideal for studying transcription factor C/EBP-delta in immune response and differentiation.

[Learn More](#)

Specification	Details
Product Name	CEBP Delta/CEBPD Rabbit pAb
Remarks/Alias	CELF; CRP3; C/EBP-delta; NF-IL6-beta; CEBP Delta/CEBPD
Species	Human
Gene ID (Human)	1052
Gene ID	1052
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human CEBP Delta/CEBPD (NP_005186.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49716
Protein Weight	28kDa
Shipping	Ice bag