

Anti-Cebp Delta/Cebpd Rabbit Polyclonal Antibody For Wb, Elisa - P49716

Item Number: CM0018571

Introduction

Rabbit polyclonal antibody targeting human CEBP Delta/CEBPD (P49716). Validated for WB and ELISA, with mouse cross-reactivity. Ideal for studying transcriptional regulation in immune and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	CEBP Delta/CEBPD Rabbit pAb
Remarks/Alias	CELF; CRP3; C/EBP-delta; NF-IL6-beta; CEBP Delta/CEBPD
Species	Human
GenelD (human)	1052
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human CEBP Delta/CEBPD (NP_005186.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P49716
Protein Weight	28kDa
Shipping	Ice bag