


## Anti-C/Ebpα Polyclonal Antibody For Wb, Elisa - P49715

Item Number: CM0023518



### Introduction

Specific rabbit polyclonal antibody against C/EBPα (CEBPA), validated for WB and ELISA. Cross-reacts with mouse. Ideal for studying myeloid differentiation, adipogenesis, and liver function.

[Learn More](#)

| Parameter        | Value                                                                                                                      |
|------------------|----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | C/EBPα Rabbit pAb                                                                                                          |
| Remarks/Alias    | CEBP; C/EBP-alpha; C/EBPα                                                                                                  |
| Species          | Human                                                                                                                      |
| GeneID (Human)   | 1050                                                                                                                       |
| GeneID           | 1050                                                                                                                       |
| Immunogen        | synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 150-250 of human C/EBPα(NP_004355.2). |
| Source           | Rabbit                                                                                                                     |
| Category         | Polyclonal Antibodies                                                                                                      |
| Application      | WB, ELISA                                                                                                                  |
| Cross Reactivity | Mouse                                                                                                                      |
| SWISS            | P49715                                                                                                                     |
| Protein Weight   | 38kDa                                                                                                                      |
| Shipping         | Ice bag                                                                                                                    |