

Anti-Cebpa Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P53566

Item Number: CM0025824

Introduction

CEBPA Rabbit polyclonal antibody from CamelBio is validated for WB, IF-P, IHC-P, ELISA. Detects CCAAT/enhancer-binding protein alpha, a transcription factor essential for myeloid differentiation, adipogenesis, and liver function.

[Learn More](#)

Parameter	Specification
Product Name	CEBPA Rabbit pAb
Remarks/Alias	Cebp; CBF-A; C/ebpalph
Species	Mouse
GeneID (Human)	1050
GeneID	12606
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-117 of mouse CEBPA (NP_031704.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53566
Protein Weight	37 kDa/36 kDa/26 kDa/41 kDa
Shipping	Ice bag