


Anti-Cebpa Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P49715

Item Number: CM0016492



Introduction

High-quality rabbit polyclonal antibody targeting human CEBPA (C/EBP-alpha). Validated for Western blot (WB), immunohistochemistry (IHC-P), and ELISA. Reacts with human, mouse, and rat. Ideal for studies on transcription factor function in myeloid differentiation, adipogenesis, and hepatogenesis.

[Learn More](#)

Parameter	Value
Product Name	CEBPA Rabbit pAb
Remarks/Alias	CEBP; C/EBP-alpha; CEBPA
Species	Human
GeneID (Human)	1050
GeneID	1050
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 111-180 of human CEBPA (NP_004355.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49715
Protein Weight	38kDa
Shipping	Ice bag