

Anti-C/Ebpb Monoclonal Antibody For Wb - P17676

Item Number: CM0003794



Introduction

C/EBPB Rabbit Monoclonal Antibody for Western blotting and ELISA. Detects human CCAAT/enhancer-binding protein beta, a transcription factor regulating immune and inflammatory responses, adipogenesis, and cell differentiation.

[Learn More](#)

Parameter	Value
Product Name	C/EBPB Rabbit mAb
Remarks/Alias	TCF5; IL6DBP; NF-IL6; C/EBP-beta; C/EBPB
Species	Human
Gene ID (Human)	1051
Gene ID	1051
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 246-345 of human C/EBPB (P17676).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P17676
Protein Weight	36kDa
Shipping	Ice bag