

Anti-C/Ebpb Polyclonal Antibody For Wb, Ihc-P, Elisa - P17676

Item Number: CM0014102



Introduction

Rabbit polyclonal antibody targeting C/EBPB (NF-IL6), validated for WB, IHC-P, ELISA. Detects human, mouse, and rat C/EBPB. Suitable for research on transcription regulation, inflammation, and adipogenesis.

[Learn More](#)

Parameter	Value
Product Name	C/EBPB Rabbit pAb
Remarks/Alias	TCF5, IL6DBP, NF-IL6, C/EBP-beta
Species	Human
GeneID (Human)	1051
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 171-216 of human C/EBPB (NP_005185.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17676
Protein Weight	36kDa
Shipping	Ice bag