

Anti-Cebpg Polyclonal Antibody For Wb - P53567

Item Number: CM0029926

Introduction

The CEBPG transcription factor regulates immune genes. This rabbit polyclonal antibody detects human and mouse CEBPG in WB, IHC-P, IF/ICC, and ELISA, enabling gene regulation studies.

[Learn More](#)

Parameter	Value
Product Name	CEBPG Rabbit pAb
Remarks/Alias	GPE1BP; IG/EBP-1; CEBPG
Species	Human
GeneID (Human)	1054
GeneID	1054
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human CEBPG (NP_001797.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P53567
Protein Weight	16kDa
Shipping	Ice bag