

Anti-Cecr5 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Bxw7

Item Number: CM0003374

Introduction

Rabbit monoclonal CECR5 antibody validated for WB, IF/ICC, and ELISA applications. Reacts with human CECR5, a widely expressed haloacid dehalogenase-like hydrolase domain-containing protein. Ideal for cat eye syndrome and hydrolase research.

[Learn More](#)

Parameter	Value
Product Name	CECR5 Rabbit mAb
Remarks/Alias	CECR5
Species	Human
Gene ID (Human)	27440
Gene ID	27440
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 324-423 of human CECR5 (Q9BXW7)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9BXW7
Protein Weight	46kDa
Shipping	Ice bag