

Anti-Cel Polyclonal Antibody For Wb And Elisa - P19835

Item Number: CM0028476



Introduction

Rabbit anti-CEL polyclonal antibody validated for WB and ELISA. Detects CEL in human, mouse, and rat samples. Suitable for studying lipid digestion and pancreatic function.

[Learn More](#)

Parameter	Value
Product Name	CEL Rabbit pAb
Remarks/Alias	BAL; FAP; BSDL; BSSL; CELL; FAPP; LIPA; CEase; MODY8; CEL
Species	Human
Gene ID (Human)	1056
Gene ID	1056
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-550 of human CEL (NP_001798.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19835
Protein Weight	79kDa
Shipping	Ice bag