

Anti-Cela2A Polyclonal Antibody For Wb, If-P, Elisa - P08217

Item Number: CM0024407

Introduction

Rabbit polyclonal antibody targeting human CELA2A (Elastase-2A). Validated for WB, IF-P, and ELISA with cross-reactivity to mouse and rat. Ideal for diabetes, insulin signaling, and glucose metabolism research.

[Learn More](#)

Parameter	Value
Product Name	CELA2A Rabbit pAb
Remarks/Alias	PE-1; AOMS4; ELA2A; CELA2A
Species	Human
GeneID (Human)	63036
GeneID	63036
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-150 of human CELA2A (NP_254275.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P08217
Protein Weight	29kDa
Shipping	Ice bag