

Anti-Celf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q92879

Item Number: CM0025945



Introduction

High-quality rabbit polyclonal antibody against human CELF1 (CUGBP1). Validated for WB, IHC-P, IP, ELISA. Cross-reacts with mouse and rat. Ideal for studying RNA-binding protein functions in alternative splicing and translation regulation.

[Learn More](#)

Parameter	Value
Product Name	CELF1 Rabbit pAb
Remarks/Alias	CUGBP; NAB50; NAPOR; CUG-BP; CUGBP1; hNab50; BRUNOL2; EDEN-BP; CELF1
Species	Human
GeneID (human)	10658
GeneID	10658
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 108-483 of human CELF1 (NP_941989.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92879
Protein Weight	52kDa
Shipping	Ice bag