

Anti-Celf5 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8N6W0

Item Number: CM0020315



Introduction

CELF5 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human and mouse samples. Targets the RNA splicing regulator BRUNOL5; useful for brain and cardiac research. – Q8N6W0

[Learn More](#)

Parameter	Value
Product Name	CELF5 Rabbit pAb
Remarks/Alias	CELF-5; BRUNOL5; BRUNOL-5; CELF5
Species	Human
Gene ID (Human)	60680
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-45 of human CELF5 (NP_001166144.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N6W0
Protein Weight	52kDa
Shipping	Ice bag