

Anti-Cenexin1/Odf2 Rabbit Monoclonal Antibody For Wb, Elisa - Q5Bjf6

Item Number: CM0011105

Introduction

Rabbit monoclonal antibody against human Cenexin1/ODF2, validated for WB and ELISA. Targets a scaffold protein in sperm tail outer dense fibers and centrosome with a predicted MW of 95 kDa. Ideal for spermatogenesis and centrosome research.

[Learn More](#)

Parameter	Value
Product Name	Cenexin1 / ODF2 Rabbit mAb
Remarks/Alias	CT134; ODF84; ODF2/1; ODF2/2; Cenexin1 / ODF2
Species	Human
GeneID (Human)	4957
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 730-829 of human Cenexin1 / ODF2 (Q5BJF6)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q5BJF6
Protein Weight	95kDa
Shipping	Ice bag