


Anti-Cenpa Rabbit Polyclonal Antibody For Wb, Elisa - P49450

Item Number: CM0020166



Introduction

Rabbit polyclonal antibody targeting human CENPA (centromere protein A), validated for WB and ELISA. Used in research on centromere identity, kinetochore assembly, and mitotic progression.

[Learn More](#)

Parameter	Value
Product Name	CENPA Rabbit pAb
Remarks/Alias	CenH3; CENP-A; CENPA
Species	Human
Gene ID (Human)	1058
Gene ID	1058
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CENPA (NP_001800.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P49450
Protein Weight	16kDa
Shipping	Ice bag