

Anti-Cenpb Rabbit Polyclonal Antibody For Wb, Elisa - P07199

Item Number: CM0021240



Introduction

Rabbit polyclonal antibody against human major centromere autoantigen B (CENPB). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for centromere and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	CENPB Rabbit pAb
Remarks/Alias	CENPB
Species	Human
GeneID (Human)	1059
GeneID	1059
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 500 to the C-terminus of human CENPB (NP_001801.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P07199
Protein Weight	65kDa
Shipping	Ice bag