


Anti-Cenpc Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q03188

Item Number: CM0012609



Introduction

High-quality rabbit monoclonal antibody targeting human CENPC (Centromere protein C, ~107 kDa). Validated for WB, IF/ICC, and ELISA applications with cross-reactivity to rat. Ideal for kinetochore and centromere research.

[Learn More](#)

Parameter	Value
Product Name	CENPC Rabbit mAb
Remarks/Alias	MIF2; hcp-4; CENP-C; CENPC1; CENPC
Species	Human
Gene ID (Human)	1060
Gene ID	1060
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 380-520 of human CENPC (Q03188)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q03188
Protein Weight	107kDa
Shipping	Ice bag