


Anti-Cenpk Polyclonal Antibody For Wb, Elisa - Q9Bs16

Item Number: CM0025405



Introduction

Anti-CENPK rabbit polyclonal antibody validated for Western blot (WB) and ELISA with human reactivity. Targets centromere protein K, a key kinetochore component essential for chromosome segregation and mitosis research.

[Learn More](#)

Parameter	Value
Product Name	CENPK Rabbit pAb
Remarks/Alias	P33; Solt; CENP-K; FKSG14; AF5alpha; CENPK
Species	Human
GeneID	64105
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BS16
Protein Weight	32kDa
Shipping	Ice bag