


Anti-Cep44 Rabbit Polyclonal Antibody For Wb, Elisa - Q9C0F1

Item Number: CM0029557



Introduction

CamelBio offers the CEP44 Rabbit Polyclonal Antibody, validated for WB and ELISA, with cross-reactivity to human and mouse CEP44. Ideal for centrosome and centriole biogenesis research. Ships on ice bag.

[Learn More](#)

Parameter	Value
Product Name	CEP44 Rabbit pAb
Remarks/Alias	PS1TP3; KIAA1712; CEP44
Species	Human
GeneID (Human)	80817
GeneID	80817
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human CEP44 (NP_001035247.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9C0F1
Protein Weight	44kDa
Shipping	Ice bag