

Anti-Cep55 Polyclonal Antibody For Wb, If/Icc, Elisa - Q53Ez4

Item Number: CM0028718



Introduction

Anti-CEP55 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA in Human, Mouse, Rat. Targets centrosomal protein of 55 kDa, key for mitotic exit and cytokinesis, ideal for cancer research.

[Learn More](#)

Parameter	Value
Product Name	CEP55 Rabbit pAb
Remarks/Alias	CT111; MARCH; URCC6; C10orf3; CEP55
Species	Human
Gene ID (Human)	55165
Gene ID	55165
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human CEP55 (NP_060601.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q53EZ4
Protein Weight	54kDa
Shipping	Ice bag