

Anti-Cep57L1 Polyclonal Antibody For Wb, Elisa - Q8Iyx8

Item Number: CM0024983

Introduction

Rabbit polyclonal anti-CEP57L1 antibody validated for WB and ELISA on human samples. Targets the centrosomal protein involved in microtubule attachment, ideal for cell division and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	Anti-CEP57L1 Polyclonal Antibody
Remarks/Alias	cep57R; C6orf182; bA487F23.2; CEP57L1
Species	Human
Gene ID (Human)	285753
Gene ID	285753
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-167 of human CEP57L1 (NP_776191.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8IYX8
Protein Weight	54kDa
Shipping	Ice bag