

Anti-Cep68 Polyclonal Antibody For Wb And Elisa - Q76N32

Item Number: CM0024597

Introduction

Rabbit polyclonal antibody targeting human CEP68, suitable for WB and ELISA. Cross-reacts with mouse and rat. Useful for investigating centrosome cohesion, CDK5RAP2 localization, and rootletin filament formation.

[Learn More](#)

Parameter	Value
Product Name	CEP68 Rabbit pAb
Remarks/Alias	KIAA0582; CEP68
Species	Human
GeneID (Human)	23177
GeneID	23177
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human CEP68 (NP_001306030.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q76N32
Protein Weight	81kDa
Shipping	Ice bag