

Anti-Cep83 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y592

Item Number: CM0030677

Introduction

High-quality rabbit polyclonal antibody against human CEP83 (CCDC41/NPHP18), validated for Western blotting and ELISA. Suitable for detecting CEP83 in human and mouse samples. Research use only.

[Learn More](#)

Parameter	Value
Product Name	CEP83 Rabbit pAb
Remarks/Alias	CCDC41; NPHP18; NY-REN-58; CEP83
Species	Human
Gene ID (Human)	51134
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-260 of human CEP83 (NP_057206.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y592
Protein Weight	83kDa
Shipping	Ice bag