

Anti-Cereblon (Crbn) Rabbit Polyclonal Antibody For Wb, Elisa - Q96Sw2

Item Number: CM0030958



Introduction

Rabbit polyclonal antibody targeting human cereblon (CRBN), validated for Western blot and ELISA. Cross-reacts with mouse and rat, suitable for studying ubiquitin-mediated degradation, neurobiology, and immunomodulation.

[Learn More](#)

Parameter	Value
Product Name	CRBN Rabbit pAb
Remarks/Alias	MRT2; MRT2A; CRBN
Species	Human
Gene ID (Human)	51185
Gene ID	51185
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human CRBN (NP_057386.2?).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96SW2
Protein Weight	51kDa
Shipping	Ice bag