

## Anti-Cerk Rabbit Polyclonal Antibody For Wb, Elisa - Q8Tct0

Item Number: CM0019109

### Introduction

Rabbit polyclonal antibody against human ceramide kinase (CERK), validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Suitable for investigating ceramide-1-phosphate signaling pathways.

[Learn More](#)

| Parameter        | Value                                                                                                                                    |
|------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CERK Rabbit pAb                                                                                                                          |
| Remarks/Alias    | LK4; hCERK; dA59H18.2; dA59H18.3; CERK                                                                                                   |
| Species          | Human                                                                                                                                    |
| Gene ID (Human)  | 64781                                                                                                                                    |
| Gene ID          | 64781                                                                                                                                    |
| Immunogen        | Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 400 to the C-terminus of human CERK (NP_073603.2). |
| Source           | Rabbit                                                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                                                    |
| Application      | WB, ELISA                                                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                        |
| SWISS            | Q8TCT0                                                                                                                                   |
| Protein Weight   | 60kDa                                                                                                                                    |
| Shipping         | Ice bag                                                                                                                                  |