

Anti-Cers5 Polyclonal Antibody For Wb And Elisa - Q8N5B7

Item Number: CM0019125

Introduction

Rabbit polyclonal antibody against human CERS5, validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for ceramide synthesis and sphingolipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	CERS5 Rabbit pAb
Remarks/Alias	Trh4; LASS5; CERS5
Species	Human
GeneID	91012
Immunogen	Synthetic peptide within amino acids 1-100 of human CERS5 (NP_001268660.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8N5B7
Protein Weight	46kDa
Shipping	Ice bag