


Anti-Ceruloplasmin Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P00450

Item Number: CM0003185



Introduction

Rabbit recombinant monoclonal antibody against human ceruloplasmin (CP). Validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for iron metabolism and copper biology studies.

[Learn More](#)

Parameter	Value
Product Name	Ceruloplasmin Rabbit mAb
Remarks/Alias	CP-2; AB073614; Ceruloplasmin
Species	Human
GeneID (Human)	1356
GeneID	1356
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 366-572 of human Ceruloplasmin (NP_000087.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00450
Protein Weight	122kDa
Shipping	Ice bag