


Anti-Ceruloplasmin Polyclonal Antibody For Wb, If/Icc, Elisa - P00450

Item Number: CM0029962



Introduction

Rabbit polyclonal antibody against human ceruloplasmin, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, rat; suitable for iron metabolism and copper transport studies.

[Learn More](#)

Parameter	Value
Product Name	Ceruloplasmin Rabbit pAb
Remarks/Alias	CP-2; AB073614; Ceruloplasmin
Species	Human
GeneID (Human)	1356
GeneID	1356
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 366-572 of human Ceruloplasmin (NP_000087.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00450
Protein Weight	122kDa
Shipping	Ice bag