


Anti-Ceruloplasmin Rabbit Polyclonal Antibody For If-P, Elisa - P00450

Item Number: CM0014311



Introduction

Rabbit polyclonal antibody targeting human Ceruloplasmin, suitable for IF-P and ELISA. Cross-reacts with mouse and rat. For studies on iron metabolism and ferroxidase activity.

[Learn More](#)

Parameter	Value
Product Name	Ceruloplasmin Rabbit pAb
Remarks/Alias	CP-2; AB073614; Ceruloplasmin
Species	Human
Gene ID (Human)	1356
Gene ID	1356
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 680-790 of human Ceruloplasmin (NP_000087.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00450
Protein Weight	122kDa
Shipping	Ice bag