


Anti-Ces1 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P23141

Item Number: CM0002085



Introduction

Rabbit monoclonal antibody against human CES1 (carboxylesterase 1). Suitable for WB, IF-P, ELISA; cross-reacts with mouse and rat. Detects 63 kDa protein. Useful for drug metabolism, lipid research, and biomarker studies.

[Learn More](#)

Parameter	Value
Product Name	CES1 Rabbit mAb
Remarks/Alias	CEH; REH; TGH; ACAT; CE-1; CES2; HMSE; SES1; HMSE1; PCE-1; hCE-1; CES1
Species	Human
Gene ID	1066
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CES1 (P23141).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23141
Protein Weight	63kDa
Shipping	Ice bag