

Anti-Ces2 Polyclonal Antibody For Wb, If/Icc, Elisa - O00748

Item Number: CM0016362

Introduction

Rabbit polyclonal antibody against human CES2 (Cocaine esterase, O00748). Validated for WB, IF/ICC, ELISA. Cross-reacts with human and mouse. Ideal for drug metabolism and detoxification studies.

[Learn More](#)

Parameter	Value
Product Name	CES2 Rabbit pAb
Remarks/Alias	iCE; CE-2; PCE-2; CES2A1; CES2
Species	Human
GeneID (Human)	8824
GeneID	8824
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 310-410 of human CES2 (NP_003860.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O00748
Protein Weight	62kDa
Shipping	Ice bag